

H1 Probe Design Documentation

Jiahao Huang

H2 Find shortest isoform for each gene in CCDS and Whole transcriptome dataset

1. Download corresponding FASTA files from NCBI FTP website
 - Human(H_sapiens) -- [CCDS](#) / [Whole transcriptome](#)
 - Mouse(M_musculus) -- [CCDS](#) / [Whole transcriptome](#)
 - Marmosets -- [GCF_000004665.1_Callithrix_jacchus-3.2_rna.fna.gz](#)
2. Find the shortest isoform for each gene
 - find_shortest_isoforms_ccds.R
 - find_shortest_isoforms_rna.R
 - find_shortest_isoforms_rna_marmosets.R (*for marmosets only*)
3. Output
 - H_sapiens_ccds_shortest_isoforms_10_21_19.fa
 - H_sapiens_rna_shortest_isoforms_10_21_19.fa
 - M_musculus_ccds_shortest_isoforms_10_21_19.fa
 - M_musculus_rna_shortest_isoforms_10_21_19.fa
 - Marmosets_rna_shortest_isoforms_10_21_19.fa

H2 PICKY 2.0

1. Input target sequence is each output FASTA file from last step (*load single file for each computation round*)
2. Parameters wt/ alternatives
 - Maximum oligo size: 46
 - Minimum oligo size: 40
 - Maximum GC content: 70
 - Minimum GC content: 30
 - Number of probes per gene: 5
 - Salt Concentration (milliM): 300
3. Parameters w/ alternatives
 1. [Maximum match length: 15 / Minimum match length: 10]
 2. [Maximum match length: 18 / Minimum match length: 15]
 3. [Maximum match length: 20 / Minimum match length: 15]
4. Output (*for each species*)

- Species_ccds_max_min.picky (ie. *H_sapiens_ccds_15_10.picky*)
- Species_rna_max_min.picky (ie. *Marmosets_rna_15_10.picky*)

noted that Marmosets does not have CCDS dataset

H2 Filtration w/ Python scripts & BBmap (Dedupe.sh)

1. probe_test_multi.py

- Input should be a folder containing all the .picky files of each species
- This script will first parse each .picky file based on the following order

CCDS -- Max: 15 Min: 10 -- #1	RNA -- Max: 15 Min: 10 -- #2
CCDS -- Max: 18 Min: 15 -- #3	RNA -- Max: 18 Min: 15 -- #4
CCDS -- Max: 20 Min: 15 -- #5	RNA -- Max: 20 Min: 15 -- #6

- Then drop all duplicated records (will keep first one)
- Remove all the records with ReverseComplement field containing continuous single nucleotide sequence longer than 5bp (ie. 'CCCCC')
- Output (for each species)
 - a log.txt file with all script running info
 - a .fa file with all probes
 - a .fq file with all probes
 - a .xlsx file with all probe records

2. dedupe.sh in=path/to/.fa out=out.fa outd=dup.fa s=1 k=20 sort=id

- This script tool from BBmap will identify duplicated probes (in this case, two probes are complementary with each other), containment probes (probe pairs where shorter one has a full length exact match with the longer one), and overlaps (in this case, minimum length of overlap is 20bp)
- In this step, every duplicated or containment probes will be stored in the dup.fa file

3. dedupe.sh in=path/to/.fa pattern=clust/cluster_%.fa fo c pto pc cc s=1 k=20 mo=20 mcs=2

- This step will generate a cluster_%.fa file for each probe set containing overlap sequence longer than 20bp

4. rm_overlap.py

- Using all the cluster_%.fa files and dup.fa file to filter the .xlsx file which contains all the probe records

- Because the probes in the original file are ordered based on their source/quality, in each overlap cluster we will keep the one with smallest index
- Output (*for each species*)
 - Species_filtered_probe.xlsx

5. head_3.py

- Pick first **3** probes for each gene in previous table
- Output (*for each species*)
 - Species_filtered_probe_head3.xlsx