

# Probe Design Documentation

---

Jiahao Huang

2023-05-09

## Find shorest isoform for each gene in CCDS and whole transcriptome dataset

---

1. Download corresponding FASTA files from [NCBI FTP](#) site
  - Human (H\_sapiens) - [whole transcriptome](#) / [CCDS](#)
  - Mouse (M\_musculus) - [whole transcriptome](#) / [CCDS](#)
  - Marmoset - [GCF\\_000004665.1\\_Callithrix\\_jacchus-3.2\\_rna.fna.gz](#) (*version used in 2019*)
2. Find the shortest isoform for each gene using the following example scripts
  - find\_shortest\_isoforms\_ccds.R
  - find\_shortest\_isoforms\_rna.R
  - find\_shortest\_isoforms\_rna\_marmosets.R (*for marmoset only*)
3. The output files of these scripts will be
  - H\_sapiens\_ccds\_shortest\_isoforms\_10\_21\_19.fa
  - H\_sapiens\_rna\_shortest\_isoforms\_10\_21\_19.fa
  - M\_musculus\_ccds\_shortest\_isoforms\_10\_21\_19.fa
  - M\_musculus\_rna\_shortest\_isoforms\_10\_21\_19.fa
  - Marmosets\_rna\_shortest\_isoforms\_10\_21\_19.fa

## Use Picky 2.0 to find probe sequences

---

1. Input target sequence is each output FASTA file from last step (*load single file for each computation round*)
2. Parameters wt/ alternatives
  - Maximum oligo size: 46
  - Minimum oligo size: 40
  - Maximum GC content: 70
  - Minimum GC content: 30
  - Number of probes per gene: 5

- `Salt Concentration (milliM): 300`

### 3. Parameters w/ alternatives

- `Maximum match length: 15 / Minimum match length: 10`
- `Maximum match length: 18 / Minimum match length: 15`
- `Maximum match length: 20 / Minimum match length: 15`

### 4. Output (for each species)

- `Species_ccds_max_min.picky` (ie. `H_sapiens_ccds_15_10.picky`)
- `Species_rna_max_min.picky` (ie. `Marmosets_rna_15_10.picky`)

\* noted that Marmosets does not have CCDS dataset

## Probe QC w/ Python scripts & BBmap (Dedupe.sh)

1. All the python functions used in this step are included in the [probe-design/scripts/2.filtration/probe.py](#) file and you can find a parsing example under the same directory titled [example.ipynb](#).
2. If you would like to run the picky parsing step on the Broad UGER cluster, related demo scripts are under [probe-design/scripts/2. filtration/cluster-run-example](#) directory.
3. Here is the original documentation in 2019 and the `probe_test_multi.py` and `rm_overlap.py` file have been moved to the `archive` folder.

#### 1. `probe_test_multi.py`

- Input should be a folder containing all the `.picky` files of each species
- The script will first parse each `.picky` file based on the following order:

<b>CCDS - Max = 15, Min = 10, #1</b>	<b>RNA - Max = 15, Min = 10, #2</b>
<b>CCDS - Max = 18, Min = 15, #3</b>	<b>RNA - Max = 18, Min = 15, #4</b>
<b>CCDS - Max = 20, Min = 15, #5</b>	<b>RNA - Max = 20, Min = 15, #6</b>

- Then drop all duplicated records (will keep first one)
- Remove all the records with ReverseComplement field containing continous single nucleotide sequence longer than 5bp (ie. 'CCCCC')
- Output (for each species)
  - a `log.txt` file with all script running info
  - a `.fa` file with all probes

- a `.fq` file with all probes
- a `.xlsx` file with all probe records

4. Once you parsed the picky files and create a dataframe with all the probes, you can filter them with the following command in [BBmap](#):

◦ `dedupe.sh in=path/to/.fa out=out.fa outd=dup.fa s=1 k=20 sort=id`

- This script tool from BBmap will identify
  - deduplicated probes - two probes are complementary with each other
  - containment probes - probe pairs where shorter one has a full length exact match with the longer one
  - overlaps - minimum length of overlap is 20bp
  - In this step, every duplicated or containment probes will be stored in the `dup.fa` file

◦ `dedupe.sh in=path/to/.fa pattern=clust/cluster_%.fa fo c pto pc cc s=1 k=20 mo=20 mcs=2`

- This step will generate a `cluster_%.fa` file for each probe set containing overlap sequence longer than 20bp

5. Then refer to the "Filtering based on bbmap (dedupe.sh)" section in [example.ipynb](#) or Run `rm_overlap.py`

- Using all the `cluster_%.fa` files and `dup.fa` file to filter the `.xlsx` file which contains all the probe records